



General Description

Junckers SYLVAstrip is designed for removal of Junckers SYLVApolish.

After treatment with Junckers SYLVAstrip, the floor should be treated with Junckers SYLVApolish.

Technical Data

Product:	Anionic tenside and nonionic tenside, soap, phosphate.
pH Value:	Concentrate: Approx 9.1 Applied Solution: Approx 8.5
Package size:	1 litre.
Coverage:	Approx. 80 - 100 m ² per litre.
Application tools:	Well wrung out cloth, or scrubbing machine.
Substrate temperature:	Minimum 5°C.
Dilution:	For 10 litre of water (25°C) the following quantities are used: Periodic cleaning: 75 ml. Basic cleaning: 250 ml
Cleaning of tools:	Water.
Storage:	Lasts for several years if unopened and stored at 20°C.

Directions for Use

Sweep, mop or vacuum clean the floor. Apply the cleaning solution as a thin even layer using a clean floor mop or a well wrung out cloth.

Rub thoroughly on very dirty areas. If necessary use a scrubbing machine. Work on small sections and rinse directly after use.

Remove the stripping solution with clean water.

If there is still some polish left, repeat the treatment.

Important!

Before use of the product, read the label on the container carefully and observe the advised precautionary safety measures.

Important!

Shake well before use.

Important!

Use as little clean water as possible to remove Junckers SYLVAstrip residue. Use a well wrung out cloth and wipes away any puddles of water on the surface immediately to prevent water damage.

Further information

Products:
F 4.8 - Junckers SYLVAstrip